

OXTM

# CUPID 6

*Portable  
Oxygen Concentrator*

Purity 87-96% @all setting

## Specification



Oxygen Flow	Pulse dose delivery system Flow settings: 1 - 6
Weight	2.0 kg ( with Slide type battery)
Dimension	175 mm (L)*70mm (W)*210 mm (H)
Noise level	39 dB(A) at setting 2
Battery	Duration more than 5 hours Charging time:150 mins
Operation	Simple control functions Easy-to-read LCD display

**Oxytek Corp.**

E-mail: [service@oxtm-o2.com](mailto:service@oxtm-o2.com)

TEL: +886 4 2321 8700

9F-2, No.239, Sec.2, Taiwan Blvd., West Dist., Taichung City, Taiwan

